

Integrated Master in Chemical Engineering

Development of formulation-specific analytical methods

Master's Thesis

By

Mafalda Sofia Pereira Brenha dos Reis

Developed in the course unit of Dissertation

Performed in

Govi Engineered Chemicals



Supervisor at FEUP: **Professor Arminda Alves**

Supervisor at GOVI: **Jan Cocquyt and Chris Van Caneyt**



Universidade do Porto
Faculdade de Engenharia

FEUP

Department of Chemical Engineering

February 2014

CONFIDENTIAL NOTICE – IMPORTANT – PLEASE READ FIRST

This document may contain confidential information proprietary to GOVI NV. It is therefore strictly forbidden to publish, cite or make public in any way this document or any part thereof without the express written permission by GOVI NV. Under no circumstance this document may be communicated to or put at the disposal of third parties; photocopying or duplicating it in any other way is strictly prohibited. Disregarding the confidential nature of this document may cause irremediable damage to GOVI NV.

All the work presented in this document is subject to a confidentiality agreement between the company, the University (FEUP) and the student. All technical information available in the company on the chemicals and methods used was made available to the author of this thesis. All the work was revised to ensure that no information pertaining to the company would be used for other purposes.

NOTA CONFIDENCIAL – IMPORTANTE – POR FAVOR LER PRIMEIRO

Este documento pode conter informação confidencial propriedade da GOVI NV. Portanto é estritamente proibido publicar, citar ou tornar público de qualquer forma este documento ou qualquer parte dele sem a expressa permissão escrita pela GOVI NV. Sob nenhuma circunstância este documento pode ser comunicado a ou posto à disposição de uma terceira parte; fotocopiar ou duplicar este de qualquer forma é estritamente proibido. Desconsiderando a natureza confidencial deste documento pode causar danos irremediáveis para a GOVI NV.

Todo o trabalho apresentado neste documento está sujeito a um acordo de confidencialidade entre a empresa, a Universidade (FEUP) e o aluno. Todas as informações técnicas disponíveis na empresa sobre os produtos químicos e os métodos utilizados foram disponibilizados para o autor desta tese. Todo o trabalho foi revisto para assegurar que nenhuma informação referente à empresa poderia ser utilizada para outros fins.

Acknowledgements

I would like to thank Jan Cocquyt for the support and help during the entire duration of the internship and for always being available to help with my thesis. I also want to thank Chris Van Caneyt for the support in the laboratory.

I would like to thank to Professor Arminda Alves for accepting my proposal for being my supervisor at FEUP and to Professor Miguel Madeira for making this internship possible.

A big thank you for my two friends Catarina Sousa and Luisa Gonçalves for being in Belgium with me and share great moments during the internship.

A special thanks to my family for being always supportive and for making this journey in college possible.

Finally, I would like thank to Tiago Marques for the support and patience, for never doubting that this is possible and for believing in me.

Abstract

Analytical methods are important to Govi since they ensure the quality of their products.

In this work five methods were tested: chloride content in emulsions, polar compounds in oils and waxes, ester and iodine number in oils and analysis of silicone oil and oil emulsion mixture. However more methods were searched and found to be possible to execute as methods that could provide us the content of phosphorus and sulphur in emulsions.

The first one was to test the chloride content in emulsion and it was proved that the best method was potentiometric with a pump to introduce AgNO_3 into the sample. The chloride content in the emulsion tested was 43.4 ± 1.2 ppm. This is the main method tested since it is the most important to Govi.

Extractable matter in the DMF and heptane phases was tested using infrared analysis after a solvent extraction. Main components found, both in waxes and oils, were unsaturated hydrocarbons, ethers and amides.

Ester and iodine numbers are used to test oils. The ester number is the difference between saponification value and acid value and both methods are based on titration. Ester number is 0.8 ± 1.0 for Oleic acid. Iodine number method is also based on titration and two different solvents were tested, CCl_4 and CHCl_3 . Results for palm oil crude sample were 56.4 ± 1.9 and $57.1 \pm 0.3 \text{ g I}_2 \cdot \text{g}^{-1}_{\text{sample}}$ for the each methods respectively.

The last method tested is very simple and provides us the content of a silicone oil and oil emulsion mixture. In this method is also used a solvent extraction. The emulsion tested has 0.53 % (m/m) of surfactants, 12.57 % (m/m) of wax and 86.91% (m/m) of silicone.

Keywords: Analytical methods, titration, infrared, solvent extraction.

Resumo

Os métodos analíticos são importantes para Govi uma vez que garantem a qualidade dos seus produtos.

Neste trabalho, foram testados cinco métodos: teor de cloreto em emulsões, compostos polares em óleos e ceras, número ésteres e iodo em óleos e a análise de uma mistura de óleo de silicone e emulsão de cera. No entanto, mais métodos foram pesquisados e conclui-se serem possíveis de executar como os métodos que podem nos fornecer o conteúdo de fósforo e enxofre em emulsões.

O primeiro método testado foi o teor de cloreto em emulsão e provou-se que o melhor método era o potenciométrico utilizando uma bomba para introduzir AgNO_3 na amostra. O teor de cloreto na emulsão foi $43,4 \pm 1,2$ ppm. Este é o teste principal método pois é o mais importante para Govi.

A matéria extraída nas fases do DMF e do heptano foi testada usando análise por infravermelhos depois de se proceder a uma extracção por solventes. Os principais componentes encontrados tanto nas ceras como nos óleos foram hidrocarbonetos insaturados, éteres e amidas.

Número de éster e de iodo são utilizados para testar óleos. O número éster é a diferença entre o valor de saponificação e valor ácido e ambos os métodos baseiam-se em titulação. O número de éster é de $0,8 \pm 1,3$ ppm para o ácido oleico. O método para testar o número de iodo é também baseada na titulação e foram testados dois solventes diferentes, CCl_4 e CHCl_3 . Os resultados para a amostra de óleo de palma bruto foram $56,4 \pm 1,9$ e $57,1 \pm 0,3 \text{ g I}_2 \cdot \text{g}^{-1}_{\text{amostra}}$. Para cada método respectivamente.

O último método testado é muito simples e fornece-nos o conteúdo da mistura de óleo de silicone e emulsão de cera. Neste método também é utilizado uma extracção por solventes. A emulsão testada tem 0,53% (m/m) de surfactantes, 12,57% (m/m) de cera e 86,91% (m/m) de silicone.

Palavras-chave: Métodos analíticos, titulação, infravermelhos, extracção por solventes.

Declaração

Declara, sob compromisso de honra, que este trabalho é original e que todas as contribuições não originais foram devidamente referenciadas com identificação da fonte.

Assinar e datar

Table of contents

| | | |
|-------|---|----|
| 1 | Introduction..... | 1 |
| 1.1 | Background and project Presentation | 1 |
| 1.2 | GOVI Engineered Chemicals..... | 1 |
| 1.3 | Contributes of the work | 1 |
| 1.4 | Thesis organization..... | 1 |
| 2 | Background and state of the art..... | 3 |
| 2.1 | Chloride content in oil mmulsion | 3 |
| 2.1.1 | Mohr’s method..... | 3 |
| 2.1.2 | Volhard’s method..... | 3 |
| 2.1.3 | Potentiometric method..... | 3 |
| 2.2 | Phosphorus content in emulsion | 4 |
| 2.3 | Sulphur content in oils and waxes | 5 |
| 2.4 | Compounds in oils and waxes extractable with a polar solvent | 6 |
| 2.5 | Ester number in bio-oils | 6 |
| 2.5.1 | Saponification value | 6 |
| 2.5.2 | Acid value | 7 |
| 2.6 | Iodine number in Bio-oils | 7 |
| 2.7 | Easy analysis method of silicone oil and wax emulsion mixture | 7 |
| 3 | Experimental work | 9 |
| 3.1 | Chloride content in Oil Emulsion..... | 9 |
| 3.1.1 | Volhard’s method..... | 9 |
| 3.1.2 | Potentiometric method..... | 10 |
| 3.2 | Compounds in oils and waxes extractable with a polar solvent..... | 11 |
| 3.3 | Ester number in bio-oils | 12 |
| 3.3.1 | Saponification value | 12 |
| 3.3.2 | Acid value | 13 |
| 3.4 | Iodine number in bio-oils | 14 |
| 3.4.1 | Carbon Tetrachloride Method | 14 |
| 3.4.2 | Chloroform Method | 15 |
| 3.5 | Easy analysis method of silicone oil and wax emulsion mixture | 16 |
| 4. | Results and Discussion | 17 |
| 4.1 | Chloride content in oil emulsion | 17 |

| | | |
|-------|---|----|
| 4.1.1 | Volhard's method | 17 |
| 4.1.2 | Potentiometric method | 18 |
| 4.2 | Compounds in oils and waxes extractable with a polar solvent | 19 |
| 4.3 | Ester number | 28 |
| 4.3.1 | Saponification value | 28 |
| 4.3.2 | Acid Value | 29 |
| 4.4 | Iodine Number | 29 |
| 4.4.1 | Carbon tetrachloride method | 29 |
| 4.4.2 | Chloroform method | 30 |
| 4.5 | Easy analysis method of silicone oil and wax emulsion mixture | 31 |
| 5 | Conclusions | 32 |
| 5.1 | Chloride content in oil emulsion | 32 |
| 5.2 | Compounds in oils and waxes extractable with a polar solvent | 32 |
| 5.3 | Ester number in bio-oils | 32 |
| 5.4 | Iodine number in bio-oils | 33 |
| 5.5 | Easy analysis method of silicone oil and wax emulsion mixture | 33 |
| 5.6 | Future Work | 33 |
| | References | 34 |
| | Appendix A – Potentiometric results | 36 |
| | Appendix B – Infrared Spectrums | 40 |
| | Appendix C – Reagents | 43 |

List of Tables

| | |
|---|----|
| Table 1: Mixtures of emulsion with 0,002 mol·L ⁻¹ NaCl and water..... | 11 |
| Table 2: Sample weights..... | 13 |
| Table 3: Guide for selecting sample size | 14 |
| Table 4: Volhard's method results on standard solutions | 17 |
| Table 5: Chloride content of emulsion components in ppm..... | 18 |
| Table 6: Extractable matter of oils A and B in both DMF and heptane phases. | 19 |
| Table 7: Extractable matter of oils A and B in both DMF and heptane phases. | 23 |
| Table 8: Saponification values of Oleic acid. | 29 |
| Table 9: Acid values of Oleic acid. | 29 |
| Table 10: Iodine values of palm oil crude with CCl ₄ method. | 30 |
| Table 11: Iodine value of palm oil crude with CHCl ₃ method. | 30 |
| Table 12: Weight of the components of wax/silicone emulsion in grams..... | 31 |
| Table 13: Concentrations of Cl ⁻ in the standards, volumes of solution titrate and of AgNO ₃ that should be spent using a pipette. | 36 |
| Table 14: Concentrations of Cl ⁻ in the standards, volumes of solution titrate and of AgNO ₃ that should be spent using a pump. | 37 |
| Table 15: Titration data and first and second derivate..... | 38 |
| Table 16: Weights of DMF, middle and heptane phases before and after drying of the oils samples. | 40 |
| Table 17: Weights of DMF, middle and heptane phases before and after drying of the oils samples. | 41 |
| Table 18: Products codes and suppliers..... | 43 |

List of Figures

| | |
|---|----|
| Figure 1: Potentiometric titration | 4 |
| Figure 2: Silver ring electrode a) Electrode used. b) Scheme of the electrode | 4 |
| Figure 3: Yellow colour proportional to phosphate concentration (From the lowest to the highest concentration). | 5 |
| Figure 4: (a) Schoniger flask ; (b) Scheme of the Schoniger flask | 6 |
| Figure 5: FTIR machine outside and inside..... | 6 |
| Figure 6: Mixture of emulsion, hexane, acetone and salt..... | 9 |
| Figure 7: a) Equipment used in the pump method; b) Capillary just above the sample surface. | 11 |
| Figure 8: Titration curves, obtained titrating with $0.01 \text{ mol}\cdot\text{L}^{-1} \text{ AgNO}_3$ sample W (green lines) and sample Y (red lines) using the pump method, potential in mV and chloride concentration in the sample in ppm..... | 18 |
| Figure 9: Infrared bonds according to wavelengths..... | 19 |
| Figure 10: IR spectrum from oil A and both separating phases, DMF and Heptane..... | 20 |
| Figure 11: Molecular structure of primary, secondary and tertiary amides..... | 20 |
| Figure 12: IF compounds analysis in the DMF phase of the separation of oil A. | 21 |
| Figure 13: IF compounds analysis in the heptane phase of the separation of oil A. | 21 |
| Figure 14: IR spectrum from oil B and both separating phases, DMF and Heptane..... | 21 |
| Figure 15: DMF, middle and heptane phases from oil B before and after drying. | 22 |
| Figure 16: IF compounds analysis in the DMF phase of the separation of oil B. | 22 |
| Figure 17: IF compounds analysis in the heptane phase of the separation of oil B. | 22 |
| Figure 18: IR spectrum from wax A and both separating phases, DMF and Heptane. | 23 |
| Figure 19: DMF, middle and heptane phases wax A before and after drying. | 23 |
| Figure 20: IF compounds analysis in the DMF phase of the separation of wax A. | 24 |
| Figure 21: IF compounds analysis in the heptane phase of the separation of wax A..... | 24 |
| Figure 22: IR spectrum from wax B and both separating phases, DMF and Heptane. | 24 |
| Figure 23: DMF, middle and heptane phases from wax B before and after drying..... | 25 |
| Figure 24: IF compounds analysis in the DMF phase of the separation of wax B..... | 25 |
| Figure 25: Ether molecular structure. | 25 |
| Figure 26: IF compounds analysis in the heptane phase of the separation of wax B..... | 25 |
| Figure 27: IR spectrum from wax C and both separating phases, DMF and Heptane. | 26 |
| Figure 28: DMF, middle and heptane phases from wax C before and after drying..... | 26 |

| | |
|--|----|
| Figure 29: IF compounds analysis in the DMF phase of the separation of wax C..... | 26 |
| Figure 30: If compounds analysis in the heptane phase of the separation of wax C. | 27 |
| Figure 31: IR spectrum from wax D and both separating phases, DMF and Heptane. | 27 |
| Figure 32: DMF, middle and heptane phases from wax D before and after drying. | 27 |
| Figure 33: IF compounds analysis in the DMF phase of the separation of wax D. | 28 |
| Figure 34: IF compounds analysis in the heptane phase of the separation of wax D..... | 28 |
| Figure 35: a) Sample before titration; b) Yellow colour appears; c) After adding starch indicator solution; d) End point of titration. | 30 |
| Figure 36: Contents in wax/silicon emulsion in %..... | 31 |
| Figure 37: Experimental titration curves (pipette method) from the standards in Table 13 relating the volume of AgNO_3 spent and the electrical potential (mV)..... | 36 |
| Figure 38: Relationship between the experimental and the theoretical Cl^- concentration in ppm in the pipette method. | 37 |
| Figure 39: Titration curves, obtained titrating with $0.01 \text{ mol}\cdot\text{L}^{-1} \text{AgNO}_3$, of sample W (green lines) and sample Y (red lines) using the weight method and chloride concentration in the sample in ppm..... | 37 |
| Figure 40: Titration curves of NaCl standard solutions (E in mV versus V AgNO_3 in ml) and experimental concentration of Cl^- in ppm. | 38 |
| Figure 41: Infrared test in Oil A. | 40 |
| Figure 42: Infrared test in oil B..... | 40 |
| Figure 43: Infrared test in wax A. | 41 |
| Figure 44: Infrared test in wax B. | 41 |
| Figure 45: Infrared test in wax C. | 42 |
| Figure 46: Infrared test in wax D..... | 42 |

Notation and Glossary

| | | |
|------------------------|--|---------------------|
| E^0 | Galvanic potential | V |
| C_{NaCl} | NaCl Concentration in standards solutions | mol·L ⁻¹ |
| $C_{TheoreticalCl-}$ | Theoretical concentration in standard solutions | ppm |
| $C_{ExperimentalCl-}$ | Experimental concentration in standard solutions | ppm |
| $m_{emulsion}$ | Weight of emulsions in sample | g |
| m_{NaCl} | Weight of NaCl in sample | g |
| m_{water} | Weight of water in sample | g |
| E | Electrical Potential | mV |
| $m_{BeforeDrying}$ | Weight of sample before drying | g |
| $m_{AfterDrying}$ | Weight of sample after drying | g |
| m_{Sample} | Weight of sample | g |
| m_{Sample} | Weight of the sample | g |
| $V_{KOHBlank}$ | Volume of KOH used in titration of the blank | ml |
| $V_{KOHSample}$ | Volume of KOH used in titration of the sample | ml |
| MW_{KOH} | Molecular weight of KOH | g·mol ⁻¹ |
| $V_{Na_2S_2O_3Blank}$ | Volume of Na ₂ S ₂ O ₃ used in the titration of the blank | ml |
| $V_{Na_2S_2O_3Sample}$ | Volume of Na ₂ S ₂ O ₃ used in the titration of the blank | ml |
| MW_I | Molecular weight of Iodine | g·mol ⁻¹ |
| N | Normalty of Na ₂ S ₂ O ₃ | mol·L ⁻¹ |

List of Chemicals

| | |
|----------------------------------|-----------------------|
| AgNO ₃ | Silver Nitrate |
| K ₂ CrO ₄ | Potassium Dichromate |
| Ag ₂ CrO ₄ | Silver Dichromate |
| KSCN | Potassium Thiocyanate |
| AgSCN | Silver Thiocyanate |
| I ₂ | Iodine |
| H ₂ O ₂ | Hydrogen Peroxide |
| HNO ₃ | Nitric Acid |
| NaCl | Sodium Chloride |
| HCl | Hydrochloric Acid |
| KOH | Potassium Hydroxide |
| CCl ₄ | Carbon Tetrachloride |
| KI | Potassium Iodide |

$\text{Na}_2\text{S}_2\text{O}_3$ Sodium Thiosulfate

CHCl_3 Chloroform

List of Acronyms

EDTA Ethylene Diamine Tetraacetic Acid

DMF Dimethylformamide

IR InfraRed

FTIR Fourier-Transform InfraRed

IPA Isopropanol

1 Introduction

1.1 Background and project presentation

Through the years, GOVI has been trying to improve their products in order to make their clients satisfied. To make sure that they produce the best products it is necessary to test them and introduce analytical methods that can provide important information to the company and the clients. The main goals of this thesis is to improve the potential of Govi's laboratory by developing methods that are reliable and specific to their products and to analyse and select raw materials to products production. It is important to determine the method that works with a certain sample and explain why other methods do not work. All the methods were implemented for the first time in the company.

The development of accurate and reliable methods that are easy to apply on water-borne process chemicals produced by Govi is an important advantage to this company.

In this work methods related to the characterization of specific chemicals and compounds will be described together with the results obtained using these methods.

1.2 GOVI Engineered Chemicals

GOVI is a Belgium company with production facilities in Belgium, Italy and Serbia. 85% of what Belgium Govi produces is exported outside Belgium. Founded in 1910, GOVI is dedicated to the manufacturing of process chemicals used in different industries, such as woodworking industry, printing, insulation products manufacturing, coatings and textile.

1.3 Contributes of the work

This work provides to GOVI the development of time-consuming methods to test their products. The development of these methods would enable the company engineers to test samples and raw materials. With the development of this thesis they will be able to just apply them on a specific product.

1.4 Thesis organization

This thesis is organized in five chapters. In the first chapter are described the objectives of this work, a brief description of the work and the company and contributes of the work.

The chapter 2 is an introduction of all methods development in the course of the internship and their fundamentals. Chapter 3 is related to the specific procedure of each method

developed in order make the methods applicable to materials used in or produced by Govi. In chapter 4 are provided the results obtained by methods and some remarks relating to the results. Chapter 5 is about the conclusions of the work as well as the future work recommended. In this chapter are described the problems faced in the work.


2 Background and state of the art

2.1 Chloride content in oil emulsion


This test was performed in order to detect if the chloride content exceeds the 20 ppm that the customer demands.

2.1.1 Mohr's method ^[1]

In the Mohr's method the solution is titrated with silver nitrate (AgNO_3) and a precipitate of silver chloride (AgCl) is formed.




As an indicator is used potassium chromate solution (K_2CrO_4) that reacts with the chloride ions which did not react with the AgNO_3 . A red-brown precipitate of silver chromate (Ag_2CrO_4) is formed when the end point is reached (reaction 2).




2.1.2 Volhard's method ^[2]


This method uses a back titration with potassium thiocyanate (KSCN) to determine the concentration of chloride ions in a solution. An excess volume of AgNO_3 is added to the solution containing chloride ions, forming a precipitate of AgCl .



Before titration, the sample is filtrated. The indicator Fe^{3+} (iron ion) is then added and the solution is titrated with potassium thiocyanate solution. Silver ions react with thiocyanate ions to form silver thiocyanate (AgSCN) precipitate.



In the beginning, the titrated is yellow but when it reacts with Fe^{3+} it forms a red complex.



The Volhard's method is used a titration and the level of chloride is determine by the end point. The end point of the titration is dark red colour.

2.1.3 Potentiometric method ^[3]

In a potentiometric method the electrochemical potential between two electrodes is measured: the indicator electrode (IE) and the reference electrode (RE). Both of them are half-cells because together they produce a certain potential.




Figure 1: Potentiometric titration ^[4].

In this method is used a combined silver ring electrode as shown in Figure 2. This is a silver/silver chloride reference electrode.





Figure 2: Silver ring electrode a) Electrode used. b) Scheme of the electrode

A AgNO_3 solution is added to the sample and the silver ions reacts with the chloride ions. An equilibrium is formed at the metal surface of the electrode that depends on the concentration of the metal ions in the sample. Metal ions are accepted by the silver ring of the electrode and silver is released into the sample solution. The equilibrium reaction is:



2.2 Phosphorus content in emulsion ^[5]

Phosphorus can occur in many forms of phosphates: orthophosphates, condensed phosphates (pyro, meta and other polyphosphates) and organic phosphates. The first step of the determination of the phosphorus content is to converse the phosphorus into orthophosphate.

Phosphorus can appear in combination with organic matter so it is necessary to choose one digestion method before determine the total phosphorus content. These methods oxidize the existent organic matter and release the orthophosphate.

There are three possible digestion methods. The perchloric acid method is recommended to difficult samples such us sediments and is a time-consuming method. The nitric acid-sulphuric acid is the simplest method and is recommended to most samples. The last possible method is the persulfate oxidation method that can be coupled with ultraviolet light to make the digestion more efficient.

The determination of the orthophosphate is done also by one of three methods. These methods are known as colorimetric method. The vanadomolybdophosphoric acid method is used for a range of 1 to 20 mg P·L⁻¹. In this method ammonium molybdate reacts with the sample to form molybdophosphoric acid and a yellow vanadomolybdophosphoric acid is formed in the presence of vanadium. The yellow colour is proportional to phosphate concentration as shown in Figure 3.



Figure 3: Yellow colour proportional to phosphate concentration (From the lowest to the highest concentration) ^[6].

The stannous chloride method is based on the formation and reduction of molybdophosphoric acid using stannous chloride. This method is more sensitive than the vanadomolybdophosphoric acid method and it is used for a range from 0.01 to 6 mg P·L⁻¹.

The ascorbic acid method ammonium molybdate and the antimony potassium tartrate react with orthophosphate and form phosphomolybdic acid.

2.3 Sulphur content in oils and waxes

Sulphur is a non-metallic element and is essential to know its content in oils and waxes because sulphur can change the behaviour on their applications. Most oils contain from 0.1 to 3 % of sulphur^[7].

Classic methods for the determination of sulphur content in many samples involve combustion of the samples and require special equipment like the Schoniger technique. Schoniger

introduced the oxygen combustion flask (Figure 4) that can be used to measure the content of several elements^[8].




Figure 4: (a) Schoniger flask ^[9] ; (b) Scheme of the Schoniger flask ^[10] .

2.4 Compounds in oils and waxes extractable with a polar solvent^[11]

Polar compounds in a sample are a good indicator of quality of the oil and the wax. To separate polar from non-polar compounds is used a solvent extraction, a method to separate compounds based on their solubility is selected such as dimethylformamide (DMF) for dissolving more polar compound and heptane for dissolving the more nonpolar compounds. After the extraction, the layers are separated as both solvents are miscible.

Samples of the two phases are dried and investigated using Fourier-Transform Infrared (FTIR) analysis. That measures the vibrations of the atoms and bonds of the molecule. Infrared (IR) radiation passes through the sample and specific wavelengths are absorbed by specific molecules present in the sample.



Figure 5: FTIR machine outside and inside.

2.5 Ester number in bio-oils^[12]

The ester number is the number of milligrams of KOH necessary for saponification of the glyceryl esters in a fat or oil and it is the difference between saponification value and acid value.

2.5.1 Saponification value^[13]

The saponification value may be defined as the number of milligrams of potassium hydroxide required to neutralize the fatty acids obtained by complete hydrolysis of one gram of oil or fat.

The expected acid value was between 196.0 and 204.0 mg KOH·g⁻¹. These values were previously determined and are part of the company's data.

2.5.2 Acid value ^[14]


The acid value may be defined as the number of milligrams of potassium hydroxide required to neutralize the free acids obtained in one gram of oil or fat.

The expected acid value was between 195.0 and 203.0 mg KOH·g⁻¹. These values were previously determined and are part of the company's data.


2.6 Iodine number in Bio-oils ^[15]

Iodine number is a measure of the unsaturation of fatty acids and is expressed in grams of iodine (I₂) per 100 grams of sample. The higher the iodine value the more double bonds are present, which consequently reflects the reactivity of the oil^[16].


The determination of this number involves 3 main steps ^[17]. The first step is the reaction of the sample with the Wijs solution in excess.



Then the excess of Wijs solution reacts with potassium iodide (KI)



and with titration with sodium thiosulphate (Na₂S₂O₃) solution it is possible to determine the amount of released iodine



The expected iodine value was between 48.0 and 58.0 g I₂·100 g⁻¹ sample. These values were previously determined and are part of the company data.

2.7 Easy analysis method of silicone oil and wax emulsion mixture

Emulsion is a dispersion of not miscible liquids, wax and silicone. They are used in many applications like cosmetics, drugs and paints.

The emulsion upon which this method was applied it is a mixture of a hydrocarbon wax emulsion and a polydimethylsiloxane (silicone) emulsion. To create an emulsion is necessary to add a surfactant (surface active agent) as a droplet stabilizing agents such as surfactants. In this analysis are used a nonionic surfactant. The goal of this method is to make a quantitative analysis of each component of the emulsion. However, it is necessary to use a solvent extraction. Heptane was used in order to separate the different components.

3 Experimental work

3.1 Chloride content in Oil Emulsion

Chloride is analysed in an emulsion. It is impossible to analyse the chloride in such a turbid and colourful sample. So it is important to use a solvent extraction and heat the emulsion to break it. Then the water phase is used to measure chloride where they are expected to be found.

3.1.1 Volhard's method^[2]

- **Reagents**

All specifications (product code and supplier) of the reagents are in Appendix C – Reagents.

1. Potassium thiocyanate
2. Silver nitrate solution ($0.1 \text{ mol}\cdot\text{L}^{-1}$)
3. Ferric ammonium sulphate ($\text{NH}_4\text{Fe}(\text{SO}_4)_2$)
4. Hydrogen peroxide 33 % (H_2O_2)
5. Nitric acid (HNO_3) concentrated

- **Procedure**

- a) **Solutions preparation**

It was necessary to prepare two solutions: KSCN solution and $\text{NH}_4\text{Fe}(\text{SO}_4)_2$ solutions. To prepare the first one it was weighted 2.43 g of solid KSCN and add to it 250 ml of distilled water in a volumetric flask. To the second one it was weighted 8 g of $\text{NH}_4\text{Fe}(\text{SO}_4)_2$ and added to it 20 ml of distilled water and a few drops of HNO_3 .

- b) **Sample preparation**

The oil emulsion was first broken by adding to 70 g of emulsion, 70 g of hexane, 70 g of acetone, and 1 g of calcium salt of formic acid.



Figure 6: Mixture of emulsion, hexane, acetone and salt.

After adding all the components to the emulsion, let it boil for about two hours. The emulsion is broken when the colour is dark brown.

After breaking the oil emulsion, it is cooled and transferred to a separating funnel.

The water phase (lowest phase) collected after breaking the emulsion is dark. So it is necessary to make it more transparent with hydrogen peroxide ^[3] to make sure that when the titration is done the colour change can be noticed.

Weight 6 g of the sample into a 500 ml flask. Add precisely 5 ml of AgNO_3 solution ($0.1 \text{ mol}\cdot\text{L}^{-1}$), 10 ml of water, 2 ml of concentrated HNO_3 and a few boiling chips. Boil the mixture in the fume hood with cooling.

Add to the solution 5 ml of 33% H_2O_2 . Keep it boiling until the colour disappears. Cool the solution, filter it and wash the solid residue with water. Transfer the sample to a 500 ml volumetric flask and fill in with water.

c) Titration

Add precisely 100 ml of sample into a conical flask and 1 ml of saturated $\text{NH}_4\text{Fe}(\text{SO}_4)_2$ as an indicator. Titrate with $0.1 \text{ mol}\cdot\text{L}^{-1}$ KSCN until dark red colour appears. Repeat the titration until obtaining 3 concordant results (titres agreeing within 0,1 ml).

3.1.2 Potentiometric method

- **Equipment**

1. Silver ring electrode (DM141-SC) from Mettler Toledo.

- **Reagents**

It is used the same AgNO_3 as in the Volhard's method. Dilutions were made to the reagent.

- **Procedure**

a) Samples preparation

Tests were performed in two samples. One is a mixture of the emulsion and water and the other a mixture of emulsion with sodium chloride (NaCl) solution so that the end emulsion contains 10 ppm more than the first one. This will provide a way to verify the reliability of the test.

Table 1: Mixtures of emulsion with $0.002 \text{ mol} \cdot \text{L}^{-1}$ NaCl and water.

| | $m_{\text{emulsion}} (\text{g})$ | $m_{\text{NaCl}} (\text{g})$ | $m_{\text{water}} (\text{g})$ |
|-----------------|----------------------------------|------------------------------|-------------------------------|
| Sample Y | 429.34 | 71.4286 | - |
| Sample W | 429.34 | - | 71.4256 |

A confidential extraction method was used to break the samples and to collect the water phase for the potentiometric method.

b) Potentiometric procedure

In order to improve the method, different ways of adding the AgNO_3 solution into the sample were used.

In all methods a well-known amount of sample between 20 to 50 ml was added from the solvent phase, to a small beaker so that the electrode goes under de liquid surface.

Pump method

In this method a pump was used to introduce $0.01 \text{ mol} \cdot \text{L}^{-1}$ AgNO_3 into the sample. In the edge of the pump was placed a syringe with a tube and a capillary placed just above the liquid surface as is shown in Figure 7.



Figure 7: a) Equipment used in the pump method; b) Capillary just above the sample surface.

3.2 Compounds in oils and waxes extractable with a polar solvent

IR is used after the sample is mixed and boiled with two solvents, heptane and DMF. More solvents were tried but their separation after boiling was not efficient because they mix after heating. This is the case of DMF and Exxsol D120, and Acetone and Exxsol D120. Solvent mixing leads to bad separation of the different compound in oil and wax. The solvents used should be immiscible both when they are cold and hot.

- **Equipment**

1. FTIR equipment (Nicolet Avatar 320 FT-IR spectrometer) from Thermo Scientific

- **Solvents**

1. Heptane
2. DMF

- **Procedure**

1. Add, approximately, 50 g of sample and 50 g of each solvent.
2. Boil them for 2 hours.
3. Cool the mixture and pour it into a separating funnel.
4. Separate the lower and upper phase. Take a middle phase to make sure the other phases are pure.
5. Dry both phases in the oven and after dry in the drying machine.
6. Test both phases in the infrared machine

3.3 Ester number in bio-oils

3.3.1 Saponification value^[18]

- **Reagents**

1. Hydrochloric acid (HCl) (5 N)
2. Potassium hydroxide in ethanol (KOH) (0.5 N)
3. Phenolphthalein
4. Ethanol 96% (v/v)

- **Procedure**

1. Weight 4 g of the sample and add 50 ml of the alcoholic KOH with a pipette.
2. Prepare a blank determination and conduct simultaneously with the sample.
3. Connect the air condensers to the flask and boil the solution gently but steadily until the sample is completely saponified. This usually requires approximately 1h for normal samples. Take care that the vapour ring in the condenser does not rise to the top of the condenser or there may be some loss.

4. After the flask and the condenser have cooled somewhat but not sufficiently to jell the contents, wash the inside of the condenser with a little distilled water. Then disconnect the flask.
5. Mix 1 g of phenolphthalein and 100 g of ethanol to make phenolphthalein indicator. Add approximately 1 ml of this indicator into the sample.
5. Titrate the solution with diluted solution of 0.5 N HCl until the pink colour has just disappeared.

3.3.2 Acid value ^[14]

The acid value may be defined as the number of milligrams of potassium hydroxide required to neutralize the free acids obtained in one gram of oil or fat.

- **Reagents**

1. Alcoholic Potassium Hydroxide 0.5 N (KOH)
2. Phenolphthalein
3. Isopropanol (IPA)

- **Procedure**

1. Weight the sample according to the following acid values expected:

Table 2: Sample weights. ^[14]

| Acid value (mg KOH·g ⁻¹) | Sample size (g) | Normality KOH (N) |
|--------------------------------------|-----------------|-------------------|
| Less than 1.0 | 4 (±0.1) | 0.1 |
| 1-25 | 2 (±0.1) | 0.1 |
| 25-50 | 1 (±0.1) | 0.1 |
| 50-150 | 3 (±0.1) | 0.5 |
| Greater than 150 | 2.5 (±0.1) | 0.5 |

Put the sample into a 125 ml Erlenmeyer flask.

2. Add 50 ml of IPA using a dispensing device and a stirring bar and mix to dissipate heating in. Allow the sample to cool.
3. Prepare a reagent blank containing 50 ml of neutralized IPA. (When using 0.5 N KOH no blank is needed).
4. Add 1g of phenolphthalein and 100 g of ethanol to make phenolphthalein indicator. Using a pipette, add 1 ml of this indicator into the sample and blank flask.
5. Titrate with KOH solution until the pink colour disappears.

3.4 Iodine number in bio-oils

3.4.1 Carbon Tetrachloride Method ^[15]

- **Reagents**

1. Carbon tetrachloride (CCl_4)
2. Wijs solution
3. Potassium iodide (KI)
4. Sodium thiosulfate solution 0,1N ($\text{Na}_2\text{S}_2\text{O}_3$)
5. Starch indicator solution (1%)

- **Procedure**

a) Solutions preparation

The potassium iodide solution should be kept in a dark bottle or prepared fresh in every attempt. To prepare 100 g of 15% KI solution weight 15 g of KI and fill with water until 100 g.

b) Titration

1. Weight the sample according to the following table, and put it into a 500 ml flask. The weight of the sample must be such that there will be an excess of Wijs solution between 100 to 150% over the amount absorbed.

Table 3: Guide for selecting sample size. ^[15]

| Iodine Value | Sample weight (g) | | Weighting Accuracy (g) |
|--------------|-------------------|-------------|------------------------|
| | 100% Excess | 150% Excess | |
| <3 | 10 | 10 | ±0.001 |
| 3 | 10.5760 | 8.4613 | 0.005 |
| 5 | 6.3460 | 5.0770 | 0.0005 |
| 10 | 3.1730 | 2.5384 | 0.0002 |
| 20 | 1.5865 | 0.8461 | 0.0002 |
| 40 | 0.7935 | 0.6346 | 0.0002 |
| 60 | 0.5288 | 0.4231 | 0.0002 |
| 80 | 0.3966 | 0.3173 | 0.0001 |
| 100 | 0.3173 | 0.2538 | 0.0001 |
| 120 | 0.2644 | 0.2115 | 0.0001 |
| 140 | 0.2266 | 0.1813 | 0.0001 |
| 160 | 0.1983 | 0.1597 | 0.0001 |
| 180 | 0.1762 | 0.1410 | 0.0001 |
| 200 | 0.1586 | 0.1269 | 0.0001 |

2. Add 20 ml of CCl_4 .
3. Pipette 25 ml of Wijs solution into the flask and shake it.
4. Prepare a blank and store the two flasks for 1 hour at $25 \pm 5^\circ\text{C}$.
5. After 1 hour add to the flasks 20 ml of KI solution and 100 ml of distilled water.
6. Titrate with $\text{Na}_2\text{S}_2\text{O}_3$ solution with constant shaking. Continue the titration until the yellow colour disappears and then add 1 ml of starch indicator solution. Titrate until the blue colour disappears.

3.4.2 Chloroform Method^[19]

- **Reagents**

1. Chloroform (CHCl_3)
2. Wijs solution
3. Potassium iodide (KI)
4. Sodium thiosulfate solution 0.1 N ($\text{Na}_2\text{S}_2\text{O}_3$)
5. Starch indicator solution (1%)

- **Procedure**

- a) **Solutions preparation**

The potassium iodide solution should be kept in a dark bottle or prepared fresh in every attempt. To prepare 100 g of 15% of KI solution weight 15 g of KI and fill with water until 100 g.

- b) **Titration**

1. Weight 0.5 g of sample.
2. Add 10 ml of CHCl_3 and warm slightly and cool for 10 minutes.
3. Pipette 25 ml of Wijs solution into the flask and shake it.
4. Prepare a blank and store the two flasks for 1 hour at $25 \pm 5^\circ\text{C}$.
5. After 1 hour add to the flasks 20 ml of KI solution and 100 ml of distilled water.
6. Titrate with $\text{Na}_2\text{S}_2\text{O}_3$ solution with constant shaking. Continue the titration until the yellow colour disappear and then add 1 ml of starch indicator solution. Titrate until the blue colour disappears.

3.5 Easy analysis method of silicone oil and wax emulsion mixture

- **Reagents**

1. Hexane

- **Procedure**

1. Before destroying the emulsion, dry it in a big beaker making a thin layer of emulsion.
2. After drying add to 30.63 g of emulsion 50 ml of hexane and mix it well.
3. The wax will not dissolve into the hexane so filter the mixture.
4. Wash the filter and the beaker with hexane to make sure that is clean.
5. Dry the wax that is in the filter.
6. Add to the filtrate 63.38 g of water to separate the silicone from the surfactants.
7. Separate the two phases with a separating funnel.
8. Dry and weight the two phases.

4. Results and Discussion

No legal limits of concentration exist for any of the methods/samples mentioned on this thesis and it was not found detection limits and precision for these methods in the matrices studied.

4.1 Chloride content in oil emulsion

4.1.1 Volhard's method

- **Calculations**

Equation 10 is the formula of calculating the chloride content in Volhard's method.

$$\frac{(V_{AgNO_3} \times C_{AgNO_3} - V_{KSCN} \times C_{KSCN} \times 5) \times MW_{Cl} \times 10^3}{10^{-3}} \quad (10)$$

- **Results**

It was only realized tests in the NaCl standards solutions.

Table 4: Volhard's method results on standard solutions

| $C_{NaCl}(mol/L)$ | $C_{Theoretical Cl^-} (ppm)$ | $C_{Experimental Cl^-} (ppm)$ |
|-------------------|------------------------------|-------------------------------|
| 0.0001 | 3.5 | -7353.5 |
| 0.0002 | 7.1 | -707.1 |
| 0.001 | 35.5 | 3114.6 |

The results of chloride content using the Volhard's method (in green) are not correct so they were excluded because it is impossible to achieve negative concentrations. This method could not be suitable for the sample used or some errors have occurred during the procedure.

4.1.2 Potentiometric method

- Results

Experimental results from chloride content are presented in Figure 8.




Figure 8: Titration curves, obtained titrating with $0.01 \text{ mol} \cdot \text{L}^{-1} \text{ AgNO}_3$ sample W (green lines) and sample Y (red lines) using the pump method, potential in mV and chloride concentration in the sample in ppm.

The chloride concentration obtained by the titration curves where calculated by the formulas on the Appendix A – Potentiometric results (Titration tables) related to the second derivate.

A correction factor of 1.4 is used. This factor takes into account the amount of water present in the original emulsion (53.78%) and the initial components amounts used in the solvent extraction.

The corrected values are $52.1 \pm 0.8 \text{ ppm}$ for sample Y and $43.4 \pm 1.2 \text{ ppm}$ for sample W. The chloride content in sample Y is $8.7 \pm 2.0 \text{ ppm}$ higher than in sample W.

Several components of the emulsion were analysed in order to determine where the chloride content came from. The same potentiometric procedure and calculations were used. In table 5 it is shown the chloride content of water, emulsifiers and gums in the emulsion. The emulsion analysed has 50% of water, 2% of emulsifier 1.1 or 1.2 and 1% of emulsifier 2.1. The emulsion (sample W) has as mentioned before has $43.4 \pm 1.2 \text{ ppm}$.

Table 5: Chloride content of emulsion components in ppm.

| Component | Chloride content in the emulsion (ppm) |
|----------------|--|
| Water | 15.3 ± 0.9 |
| Emulsifier 1.1 | 16.5 ± 1.0 |
| Emulsifier 1.2 | 22.5 ± 0.2 |
| Emulsifier 2.1 | 2.0 ± 0.2 |

The emulsion tested has 43.4 ± 1.2 ppm. According to this test, if the emulsifier 1.1 was used the emulsion has 33.8 ± 2.1 ppm and if the emulsifier 1.2 was used than the emulsion should have 39.8 ± 1.3 ppm. The components mentioned in this test are just 53% of the emulsion. The missing content of chloride comes probably from other components.

4.2 Compounds in oils and waxes extractable with a polar solvent

• Results

This test was carried out in several samples of oil and samples.

The infrared spectrum is from 4000 to 400 cm^{-1} . This is a mid-infrared region that can be divided into four regions.

However, it is possible to analyse the spectrum, sometimes it is impossible to be sure that a certain bond exists or not. Certain peaks are difficult to analyse because they are not exactly as theory of infrared says.




Figure 9: Infrared bonds according to wavenumbers. ^[20]

Oils

The results of extractable matter of both oils A and B are presented in Table 6.

Table 6: Extractable matter of oils A and B in both DMF and heptane phases.

| | Extractable matter (%) | |
|-------|------------------------|---------------|
| | DMF Phase | Heptane Phase |
| Oil A | 4.87 | 95.13 |
| Oil B | 0.74 | 99.26 |

Sample: Oil A

An infrared analysis was performed in oil A and in both extractable phases as is possible to see in Figure 10.




Figure 10: IR spectrum from oil A and both separating phases, DMF and Heptane.

In the DMF phase of the separation (Figure 12) it is possible to identify hydrocarbons because the stretch C-H bond is visible between $3000\text{--}2850\text{ cm}^{-1}$. It is a saturated hydrocarbon because the C-H wavelength is lower than 3000 cm^{-1} . It is also visible the bending C-H bond around $1450\text{--}1375\text{ cm}^{-1}$. The equipment also identifies a tertiary amides existence but the amide existent is not tertiary because this kind of amide is known for not having an N-H bond as is possible to see in the spectrum. However, this bond is visible between $3500\text{--}3100\text{ cm}^{-1}$. The C=O bond is visible around 1670 cm^{-1} .




Figure 11: Molecular structure of primary, secondary and tertiary amides.

In Figure 12, C=C bond is visible around 1600 cm^{-1} . This bond and the C-H bond found around $3000\text{--}3100\text{ cm}^{-1}$ help to identify an alkene.




Figure 12: IR compounds analysis in the DMF phase of the separation of oil A.

In the heptane phase of the separation, it is also possible to see a saturated hydrocarbon. The bending and stretch C-H bonds are visible on the same wavelength of the lower phase. However, as it is possible to see in Figure 13, FTIR analysis is not correct because it identifies carboxylic acids existence. For it exist, a C=O bond around $1725\text{--}1700\text{ cm}^{-1}$ should be visible.




Figure 13: IR compounds analysis in the heptane phase of the separation of oil A.

Sample: Oil B

An infrared analysis was performed in oil B and in both extractable phases as is possible to see in Figure 14.




Figure 14: IR spectrum from oil B and both separating phases, DMF and Heptane.



Figure 15: DMF, middle and heptane phases from oil B before and after drying.

In Figure 16 it is identified an amide but the C=O bond should be more visible between $1670\text{--}1640\text{ cm}^{-1}$. So it is not possible to find an amide in this phase. It is possible to identify an alkene because a C=C above 1600 cm^{-1} . It also exists a stretch and bending C-H bonds between $3000\text{--}2850\text{ cm}^{-1}$ and $1450\text{--}1375\text{ cm}^{-1}$, respectively, characteristics of a saturated hydrocarbon.




Figure 16: IR compounds analysis in the DMF phase of the separation of oil B.

This is a simple spectrum analysis. In Figure 17 it is possible to identify a saturated hydrocarbon with a stretch and bending C-H bond between $3000\text{--}2850\text{ cm}^{-1}$ and $1450\text{--}1375\text{ cm}^{-1}$, respectively.




Figure 17: IR compounds analysis in the heptane phase of the separation of oil B.

Waxes

The results of extractable matter of waxes A, B, C and D are presented in Table 7.

Table 7: Extractable matter of oils A and B in both DMF and heptane phases.

| | Extractable matter (%) | |
|--------------|------------------------|---------------|
| | DMF Phase | Heptane Phase |
| Wax A | 0.46 | 99.54 |
| Wax B | 1.86 | 98.14 |
| Wax C | 0.94 | 99.06 |
| Wax D | 0.38 | 99.02 |

Sample: Wax A

An infrared analysis was performed in wax A and in both extractable phases, as it is possible to see in Figure 18.




Figure 18: IR spectrum from wax A and both separating phases, DMF and Heptane.



Figure 19: DMF, middle and heptane phases wax A before and after drying.

In Figure 20 it is also visible a saturated hydrocarbon with a stretch and bending C-H bonds. An error occurs in the FTIR analysis because sulfoxides bonds are not visible. The S=O bond should be visible around 1050 cm^{-1} and should be a stronger and more visible bond. However, this spectrum has more picks than the normal hydrocarbon spectrums. This could mean that more compounds exist.




Figure 20: IF compounds analysis in the DMF phase of the separation of wax A.

Figure 21 is also very easy to analyse because it is a characteristic saturated hydrocarbon bond. The bonds found are the same as the ones found in the others saturated hydrocarbon spectrums.




Figure 21: IF compounds analysis in the heptane phase of the separation of wax A.

Sample: Wax B

An infrared analysis was performed in wax B and in both extractable phases as is possible to see in Figure 22.




Figure 22: IR spectrum from wax B and both separating phases, DMF and Heptane.



Figure 23: DMF, middle and heptane phases from wax B before and after drying.

In Figure 24 it is also possible to identify a saturated hydrocarbon. This is a long chain hydrocarbon because it is visible a peak around 720 cm^{-1} . However, this spectrum has some peaks impossible to identify and this peaks could mean that other compounds exist.




Figure 24: IF compounds analysis in the DMF phase of the separation of wax B.

In Figure 26 it is found a C-O bond around 1115 cm^{-1} that characterizes ether. As in most spectrums an unsaturated hydrocarbon is also identified.




Figure 25: Ether molecular structure.




Figure 26: IF compounds analysis in the heptane phase of the separation of wax B.

Sample: Wax C

An infrared analysis was performed in wax C and in both extractable phases as is possible to see in Figure 27.




Figure 27: IR spectrum from wax C and both separating phases, DMF and Heptane.



Figure 28: DMF, middle and heptane phases from wax C before and after drying.

In Figure 29 it is visible a C=O bond around 1715 cm^{-1} that characterizes a ketone. In this figure and in Figure 30 it is visible a long chain saturated hydrocarbon.




Figure 29: IF compounds analysis in the DMF phase of the separation of wax C.




Figure 30: IR compounds analysis in the heptane phase of the separation of wax C.

Sample: Wax D

An infrared analysis was performed in wax D and in both extractable phases as is possible to see in Figure 31.




Figure 31: IR spectrum from wax D and both separating phases, DMF and Heptane.



Figure 32: DMF, middle and heptane phases from wax D before and after drying.

In the DMF phase spectrum it is visible also a long chain hydrocarbon.




Figure 33: IF compounds analysis in the DMF phase of the separation of wax D.

In this spectrum (Figure 34) it is found a C-O bond around 1115 cm^{-1} that characterizes ether. As in most spectrums a saturated hydrocarbon is also identified.




Figure 34: IF compounds analysis in the heptane phase of the separation of wax D.

4.3 Ester number

4.3.1 Saponification value

- Calculations

$$\text{Saponification value} = \frac{MW_{\text{KOH}} \times N_{\text{HCl}} \times (V_{\text{HCl Blank}} - V_{\text{HCl Sample}})}{m_{\text{sample}}} \quad (11)$$

$$MW_{\text{KOH}} = 56.1\text{ g}\cdot\text{mol}^{-1}$$

- Results

Sample: Oleic acid

The standardization of HCl was made with 0.5 N KOH. The concentration of HCl is $0.4717\text{ mol}\cdot\text{L}^{-1}$.

Table 8: Saponification values of Oleic acid.

| $m_{\text{sample}} \text{ (g)}$ | $V_{\text{KOH Blank}} \text{ (ml)}$ | $V_{\text{KOH Sample}} \text{ (ml)}$ | Saponification value (mg KOH·g⁻¹) |
|-------------------------------------|-------------------------------------|--------------------------------------|---|
| 4.05 | 53.0 | 21.8 | 203.9 |
| 4.09 | 53.0 | 21.7 | 202.5 |
| 4.07 | 53.0 | 21.9 | 202.2 |
| Average ± Standard Deviation | | | 202.9 ± 0.9 |

4.3.2 Acid Value

- **Calculations**

$$\text{Acid value} = \frac{(V_{\text{KOH sample}} - V_{\text{KOH blank}}) \times N_{\text{KOH}} \times MW_{\text{KOH}}}{m_{\text{sample}}} \quad (12)$$

Using 0.5 N KOH

$$\text{Acid value} = \frac{V_{\text{KOH sample}} \times N_{\text{KOH}} \times MW_{\text{KOH}}}{m_{\text{sample}}} \quad (13)$$

$MW_{\text{KOH}} = 56.1 \text{ g} \cdot \text{mol}^{-1}$

- **Results**

Sample: Oleic acid

Table 9: Acid values of Oleic acid.

| $m_{\text{sample}} \text{ (g)}$ | $V_{\text{KOH Sample}} \text{ (ml)}$ | Acid value (mg KOH·g⁻¹) |
|-------------------------------------|--------------------------------------|---|
| 2.50 | 21.8 | 202.0 |
| 2.52 | 18.2 | 202.6 |
| 2.53 | 18.2 | 201.8 |
| Average ± Standard deviation | | 202.1±0.4 |

4.4 Iodine Number

4.4.1 Carbon tetrachloride method

- **Calculation**

$$\text{Iodine value} = \frac{(V_{\text{Na}_2\text{S}_2\text{O}_3 \text{ Blank}} - V_{\text{Na}_2\text{S}_2\text{O}_3 \text{ Sample}}) \times N_{\text{Na}_2\text{S}_2\text{O}_3} \times MW_{\text{I}}}{m_{\text{sample}}} \times \frac{100}{1000} \quad (14)$$

$MW_{\text{I}} = 126.9 \text{ g} \cdot \text{mol}^{-1}$

- Results




Figure 35: a) Sample before titration; b) Yellow colour appears; c) After adding starch indicator solution; d) End point of titration.

Sample: Palm Oil crude

The experimental average iodine number is $56.4 \text{ g I}_2 \cdot 100 \text{ g}^{-1} \text{ sample}$.

Table 10: Iodine values of palm oil crude with CCl_4 method.

| $V_{\text{Na}_2\text{S}_2\text{O}_3 \text{ Blank}} (\text{ml})$ | $V_{\text{Na}_2\text{S}_2\text{O}_3 \text{ Sample}} (\text{ml})$ | $m_{\text{sample}} (\text{g})$ | Iodine value ($\text{g I}_2 \cdot 100 \text{g}^{-1} \text{ sample}$) |
|---|--|--------------------------------|--|
| 48.9 | 18.5 | 0.70 | 55.1 |
| 49.1 | 17.7 | 0.69 | 57.7 |
| Average \pm Standard Deviation | | | 56.4 \pm 1.9 |

4.4.2 Chloroform method

- Calculation^[16]

$$\text{Iodine value} = \frac{(V_{\text{Na}_2\text{S}_2\text{O}_3 \text{ Blank}} - V_{\text{Na}_2\text{S}_2\text{O}_3 \text{ Sample}}) \times N_{\text{Na}_2\text{S}_2\text{O}_3} \times \text{MW}_I}{m_{\text{sample}}} \times \frac{100}{1000} \quad (15)$$

- Results

Sample: Palm Oil crude

The experimental average iodine number is $57.4 \text{ g I}_2 \cdot 100 \text{ g}^{-1} \text{ sample}$.

Table 11: Iodine value of palm oil crude with CHCl_3 method.

| $V_{\text{Na}_2\text{S}_2\text{O}_3 \text{ Blank}} (\text{ml})$ | $V_{\text{Na}_2\text{S}_2\text{O}_3 \text{ Sample}} (\text{ml})$ | $m_{\text{sample}} (\text{g})$ | Iodine value ($\text{g I}_2 \cdot 100 \text{g}^{-1}$) |
|---|--|--------------------------------|---|
| 0.51 | 29.7 | 52.9 | 57.7 |
| 0.52 | 29.4 | 52.9 | 57.3 |
| 0.52 | 29.5 | 52.9 | 57.1 |
| Average \pm Standard Deviation | | | 57.4 \pm 0.3 |

4.5 Easy analysis method of silicone oil and wax emulsion mixture

Only one test was performed using this method so no standard deviation is calculated. The results are in mass percentage.

Table 12: Weight of the components of wax/silicone emulsion in grams.

| | Surfactants | Wax | Silicone |
|---------------------------|-------------|------|----------|
| <i>m</i> after drying (g) | 0.15 | 3.62 | 25.03 |




Figure 36: Contents in wax/silicon emulsion in %.

5 Conclusions

All methods tested and implemented were effective and allowed to characterize and test Govi's products and raw materials.

5.1 Chloride content in oil emulsion

Different methods were tested to determine chloride content in the oil emulsion selected. However, only two methods provided good results.

Volhard's method lacked accuracy and was therefore not suitable for these samples.

The potentiometric method using a silver ring electrode and the syringe pump to transfer the AgNO_3 to the samples was more sensitive and can be applied to coloured even turbid samples.

In one emulsion sample 10 ppm Cl^- was added. Then the original sample was compared to this sample with higher chloride concentration, using the latter method, 8.7 ppm was traced back, showing that this method is reliable.

5.2 Compounds in oils and waxes extractable with a polar solvent

Both polar and nonpolar solvents were selected and added to the wax or oil. After boiling, polar compounds are concentrated in the polar phase (DMF phase) spectrum and nonpolar compounds on the nonpolar phase (heptane phase). Infrared analysis provides mainly qualitative information of chemical bonds.

Quantitative information on polar compounds can be obtained by weighting the amounts in each solvent. With this method it is possible to determine the extractable matter in the two solvent phases. Wax B showed a higher amount of DMF extractable matter than the others.

5.3 Ester number in bio-oils

The ester number is the difference between saponification and acid values. The average saponification value is $202.9 \pm 0.88 \text{ mg KOH} \cdot \text{g}^{-1}$ and the average acid value is $202.1 \pm 0.42 \text{ mg KOH} \cdot \text{g}^{-1}$. So the ester number is $0.8 \pm 1.0 \text{ mg KOH} \cdot \text{g}^{-1}$.

5.4 Iodine number in bio-oils

For the palm oil crude sample the iodine value in the technical data sheet was between 48.0 and 58.0 g I₂·100 g⁻¹_{sample} whereas on the work was found 56.4±1.87 and 57.1±0.31 g I₂·100 g⁻¹_{sample} for CCl₄ and CHCl₃ methods respectively.

5.5 Easy analysis method of silicone oil and wax emulsion mixture

The emulsion tested has 0.52 % (m/m) of surfactants, 12.57 % (m/m) of wax and 86.91 % (m/m) of silicone. However this result might not be conclusive because only one test was performed due to the lack of time.

5.6 Future Work

In the future, the phosphorus and sulphur determination methods should be applied because they were analysed and their execution is possible. Another possible future work is to apply the extracted wax in boards in order to see if the performance improved.

One important task is to continue testing emulsion components and test new components in order to achieve the 20 ppm demanded by the costumers.

References

1. University of Canterbury. *Determination of Chloride Ion Concentration by Titration (Mohr's Method)*. www.outreach.canterbury.ac.nz. Retrieved September 22, 2013, from: http://www.outreach.canterbury.ac.nz/chemistry/documents/chloride_mohr.pdf.
2. University of Canterbury. *Determination of Chloride Ion Concentration by Titration (Volhard's Method)*. www.outreach.canterbury.ac.nz. Retrieved September 22, 2013, from: http://www.outreach.canterbury.ac.nz/chemistry/documents/chloride_volhard.pdf.
3. Metrohm. *Basics of Potentiometry*. www.metrohm.com. Retrieved October 17, 2013, from: http://www.metrohm.com/com/downloads/Support_pdf/BASICS_OF_POTENTIOMETRY.pdf.
4. *Potentiometric Titration Apparatus*. www.clipartist.net. Retrieved December 18, 2013, from: <http://clipartist.net/svg/potentiometric-titration-apparatus-september-2011-clip-art-svg-opensclipart-org-commons-wikimedia-org-clipartist-net/>.
5. Clescer, Lenore S.. *Standard Methods for the Examination of Water and Wastewater*. 20th Edition. American Public Health Association. 1998.
6. Dublin City University. www.dcu.ie. Retrieved October 15, 2013, from: <http://www.dcu.ie/chemistry/asg/mcgrawc/index.shtml>.
7. Argonne National Laboratory. *Sulphur in heavy oils, oil sands and oil shales*. web.anl.gov. Retrieved January 15, 2014 from: http://web.anl.gov/PCS/acsfuel/preprint%20archive/Files/22_3_MONTREAL_06-77_0086.pdf.
8. Levaggi, D.A., Feldstein, M.. *A rapid method for the determination of sulphur in fuel oil by the Schoniger oxygen flask method*. Journal of air pollution control association. 2012.
9. Virtuellen Museum der Wissenschaft. *Schoeniger kolben*. www.amuseum.de. Retrieved in October 20, 2013, from: <http://www.amuseum.de/chemie/Schoeniger/schoenigerFlask.htm>.
10. U.S. Pharmacopeia. *Oxygen flask combustion*. www.pharmacopeia.cn. Retrieved October 20, 2013, from: http://www.pharmacopeia.cn/v29240/usp29nf24s0_c471.html.
11. Stuart, B. (2004). *Infrared spectroscopy: fundamentals and applications*. John Wiley & Sons, Ltd. p. 1-18.
12. Lookchem. *Ester number*. www.lookchem.com. Retrieved October 17, 2013, from: http://www.lookchem.com/Dictionary/Ester_Number/.

13. Virtual Amrita laboratories universalizing education. *Estimation of saponification value for fats/oils*. amrita.vlab.co.in. Retrieved in October 10, from: <http://amrita.vlab.co.in/?sub=3&brch=63&sim=688&cnt=1>.
14. Lubrizol. *Acid Value*. www.lubrizol.com. Retrieved October 20, 2013, from: http://www.lubrizol.com/Personal-Care/Documents/Test-Procedures/Emulsifiers/TP-AATM_109-01.pdf.
15. ASTM Standard D 5554-95, *Standard Test Method for Determination of the Iodine Value of Fats and Oils*, ASTM International, West Conshohocken, PA 19428-2959.
16. Seatons. *Iodine Value*. www.seatons-uk.co.uk. Retrieved October 20, 2013, from: <http://www.seatons-uk.co.uk/home.aspx?s=62&r=109&p=507>
17. QuadraChem Laboratories. *Iodine Value of Animal and Vegetables Fats and Oils*. www.qclscientific.com. Retrieved November 7, 2013, from: <http://www.qclscientific.com/pdfs/Apps%20Notes/Iodine%20Value%20of%20Animal%20and%20Vegetable%20Fats%20and%20Oils.pdf>.
18. ASTM Standard D 5558-95, *Standard Test Method for Determination of the Saponification Value of Fats and Oils*, ASTM International, West Conshohocken, PA 19428-2959.
19. Classle. *Estimation of Iodine Value of Oil*. www.classle.net. Retrieved January 5, 2014, from: <https://www.classle.net/book/estimation-iodine-value-oil#>.
20. Stuart, B. (2004). *Infrared spectroscopy: fundamentals and applications*. John Wiley & Sons, Ltd. p. 115.

Appendix A – Potentiometric results

a) Pipette method / Weight method

The potentiometric method was tested using also a pipette to transfer AgNO_3 to the sample and also weighting AgNO_3 and transfer into the sample.

Table 13: Concentrations of Cl^- in the standards, volumes of solution titrate and of AgNO_3 that should be spent using a pipette.

| Standards | $C_{\text{Cl}^-}(\text{ppm})$ | $C_{\text{AgNO}_3}(\text{mol} \cdot \text{L}^{-1})$ | $V_{\text{NaCl solution titrated}}(\text{ml})$ | $V_{\text{AgNO}_3 \text{ Expected}}(\text{ml})$ |
|-----------|-------------------------------|---|--|---|
| 1 | 70.9 | 0.001 | 20 | 40 |
| 2 | 35.5 | 0.001 | 50 | 50 |
| 3 | 17.4 | 0.001 | 40 | 20 |
| 4 | 7.1 | 0.001 | 30 | 6 |

The use of the second derivate to determine the equivalence point is impossible because the graphic has a lot of noise. However it is possible to see in Figure 37 that the expected values of AgNO_3 volumes are achieved.




Figure 37: Experimental titration curves (pipette method) from the standards in Table 13 relating the volume of AgNO_3 spent and the electrical potential (mV).




Figure 38: Relationship between the experimental and the theoretical Cl^- concentration in ppm in the pipette method.




Figure 39: Titration curves, obtained titrating with $0.01 \text{ mol} \cdot \text{L}^{-1} \text{AgNO}_3$, of sample W (green lines) and sample Y (red lines) using the weight method and chloride concentration in the sample in ppm.

b) Pump method

Standards solutions of NaCl were prepared by dilution of $0.1 \text{ mol} \cdot \text{L}^{-1}$ and tested using a potentiometric method adding AgNO_3 with a pump.

Table 14: Concentrations of Cl^- in the standards, volumes of solution titrate and of AgNO_3 that should be spent using a pump.

| Standards | $C_{\text{Cl}^-} (\text{ppm})$ | $C_{\text{AgNO}_3} (\text{mol} \cdot \text{L}^{-1})$ | $V_{\text{NaCl solution titrated}} (\text{ml})$ | $V_{\text{AgNO}_3 \text{ Expected}} (\text{ml})$ |
|-----------|--------------------------------|--|---|--|
| 1 | 70.90 | 0.001 | 20 | 40 |

| | | | | |
|---|-------|-------|----|----|
| 2 | 35.45 | 0.001 | 20 | 20 |
| 3 | 17.43 | 0.001 | 20 | 10 |
| 4 | 7.09 | 0.001 | 30 | 6 |




Figure 40: Titration curves of NaCl standard solutions (E in mV versus $V \text{ AgNO}_3$ in ml) and experimental concentration of Cl^- in ppm.

c) Titration tables

The potential E and the volume of AgNO_3 (V) are the data achieved in the titration. $\Delta E/\Delta V$ (Equation 16) and $\Delta^2 E/\Delta V^2$ (Equation 17) are the first and second derivate. With the second derivate is possible to determine the $V \text{ AgNO}_3$ used. Table 15 it is an example of a table used in all potentiometric procedures. This table has the data of the standard with 35.45 ppm of Cl^- .

$$\frac{\Delta E}{\Delta V} = \frac{E_2 - E_1}{V_2 - V_1} \quad (16)$$

$$\frac{\Delta^2 E}{\Delta V^2} = \frac{\left(\frac{\Delta E}{\Delta V}\right)_2 - \left(\frac{\Delta E}{\Delta V}\right)_1}{V_2 - V_1} \quad (17)$$

Table 15: Titration data and first and second derivate.

| $E \text{ (mV)}$ | $V_{\text{AgNO}_3} \text{ (ml)}$ | $\Delta E/\Delta V$ | $\Delta^2 E/\Delta V^2$ |
|------------------|----------------------------------|---------------------|-------------------------|
| 11 | 0 | | |
| | | 1 | |
| 16 | 5 | | 0 |
| | | 1 | |
| 21 | 10 | | 0.2 |
| | | 2 | |

| | | |
|-----|-----|-------|
| 31 | 15 | 1.25 |
| | 4.5 | |
| 40 | 17 | 2.5 |
| | 7 | |
| 47 | 18 | 1 |
| | 8 | |
| 55 | 19 | 6 |
| | 14 | |
| 69 | 20 | 10 |
| | 24 | |
| 93 | 21 | -4 |
| | 20 | |
| 113 | 22 | -3 |
| | 17 | |
| 130 | 23 | -5.5 |
| | 6 | |
| 142 | 25 | -0.56 |
| | 3.2 | |
| 158 | 30 | -0.16 |
| | 2.4 | |
| 170 | 35 | -0.24 |
| | 1.2 | |
| 176 | 40 | |

To calculate the volume of AgNO_3 used to in the equivalence point it is necessary to calculate the point where the second derivative is zero. So an interpolation between the values 10 e -4 was done. With the volume of AgNO_3 , the volume of initial sample and the concentration of AgNO_3 it is possible to determine the concentration of Cl^- in the sample.

Appendix B – Infrared Spectrums

Table 16: Weights of DMF, middle and heptane phases before and after drying of the oils samples.

| Oil | Phase | $m_{\text{Before drying}} (g)$ | $m_{\text{After drying}} (g)$ |
|-----|---------|--------------------------------|-------------------------------|
| A | DMF | 46.50 | 2.40 |
| | Middle | 3.43 | 1.14 |
| | Heptane | 94.74 | 47.09 |
| B | DMF | 45.14 | 1.31 |
| | Middle | 1.34 | 0.44 |
| | Heptane | 102.91 | 51.11 |

Sample: Oil A




Figure 41: Infrared test in Oil A.

Sample: Oil B






Figure 42: Infrared test in oil B.

Table 17: Weights of DMF, middle and heptane phases before and after drying of the oils samples.

| Wax | Phase | $m_{\text{Before drying (g)}}$ | $m_{\text{After drying (g)}}$ |
|-----|---------|--------------------------------|-------------------------------|
| A | DMF | 41.83 | 0.20 |
| | Middle | 3.60 | 0.62 |
| | Heptane | 94.31 | 44.62 |
| B | DMF | 42.68 | 0.58 |
| | Middle | 1.68 | 0.52 |
| | Heptane | 91.75 | 46.34 |
| C | DMF | 49.59 | 0.95 |
| | Middle | 1.72 | 0.81 |
| | Heptane | 97.82 | 49.33 |
| D | DMF | 42.11 | 0.50 |
| | Middle | 1.64 | 0.52 |
| | Heptane | 108.43 | 53.17 |
| E | DMF | 40.20 | 0.20 |
| | Middle | 1.37 | 0.17 |
| | Heptane | 104.47 | 53.48 |

Sample: Wax A**Figure 43:** Infrared test in wax A.**Sample: Wax B****Figure 44:** Infrared test in wax B.

Sample: Wax C




Figure 45: Infrared test in wax C.

Sample: Wax D




Figure 46: Infrared test in wax D.

Appendix C – Reagents

Table 18: Products codes and suppliers.

| Reagent | Concentration | Code | Supplier |
|---|-------------------------|--------------|------------|
| KSCN | - | 27035.230 | Prolabo |
| AgNO ₃ | 0.1 mol·L ⁻¹ | 1.09081.1000 | Merck |
| NH ₄ Fe(SO ₄) ₂ | - | 1.03776.0500 | Merck |
| H ₂ O ₂ | 33% | 23613.297 | Merck |
| HNO ₃ | Concentrated | 9606-34 | J.T. Baker |
| NaCl | 0.1 mol·L ⁻¹ | 35616 | Alfa Aesar |
| Heptane | - | 24548.367 | Prolabo |
| DMF | - | 0464-500ML | Ambresco |
| HCl | 5 N | 30018.320 | Prolabo |
| KOH in ethanol | 0.5 N | 31308.297 | Prolabo |
| Phenolphthalein | - | 1.07233.0100 | Merck |
| Ethanol | 96% (v/v) | 84836.360 | Prolabo |
| CCl ₄ | - | 22521.293 | Prolabo |
| Wijs solution | - | 2812.1000 | Prolabo |
| KI | | 1.05043.1000 | Merck |
| Na ₂ S ₂ O ₃ | 0.1 N | 1.09147.1000 | Merck |
| Starch indicator | 1% | 1085103 | Reagecon |
| CHCl ₃ | - | 22706.326 | Prolabo |
| Hexane | - | 24611.366 | Prolabo |